




Easydrip™ Classic PEN NEEDLES

Easydrip™ Classic can be used with an injection pen for subcutaneous injection of insulin, GLP-1RA , growth hormone, etc.

Common GLP-1RA (Glucagon-Like Peptide-1 Receptor Agonists) drugs include: Semaglutide, Liraglutide, Dulaglutide, Lixisenatide, Exenatide, Albiglutide.

Ordering Information

Reorder #	Color	Needle Gauge	Needle Length	Standard Package
01-3204		32G needle	4mm	100 Pcs/Box 2,000 Pcs/Case
01-3105		31G needle	5mm	
01-3106		31G needle	6mm	

1. Needle gauge refers to the size of the needle. The higher the gauge, the thinner the needle;
2. Please contact your sales person if you have special requests for an unique size.



MDSAP



MDL

510K

CE

FMR

Austenite Extra-Thin-Walled needle increases flow and shortens injection time.




● Safety and comfort

Medical silicone lubrication of both ends of the needle improves safety and comfort.




The needle tip pierces the skin smoothly



The needle tail pierces the rubber stopper without coring or production of tiny particles.

● Painless experience

A resilient and ultra-sharp needle reduces pain and increases injection comfort.¹



- Cobalt-60 radiation sterilization leaves no chemical residue and broad-spectrum sterilization, ensuring a clean and safe injection.

STERILE R Clean & Green
Sterilized using irradiation

Reference:

1. Study on the influence of structure design and usage method on the performance of insulin pen needles